


GOOD HEALTH STARTS WITH **GOOD TEETH**

The ultrasonic toothbrush for your pet



emmi[®] PET

Did you know

that dental problems are quite common among cats and dogs? From the age of three years, eight out of ten pets have problems regarding their teeth or gums.



1.

Good health begins in the mouth



This is not only true for us humans, but also for cats and dogs. Clean, fresh breath and shining white teeth are good for animals, too. If teeth and gums are thoroughly cared for, our four-legged friends can chew properly and enjoy their food. However, if something is wrong in their mouth,

the owners often do not notice it: Our precious pets do not stop eating and they also cannot tell us, their owners, that something is wrong. In the end a professional treatment of dental calculus, caries and resulting complications is necessary – off to the animal dentist.


Ideal for health teeth and gums

What do we do to care for our teeth? We brush them! Brushing your teeth regularly is the most effective way to keep your mouth and teeth in good health – the same applies to our four-legged friends. Even though we have special animal feed for dental cleaning, for example special chewing bones, leftovers from food and bacteria still remain in the mouth. Thorough cleaning and preventative care are done best with toothbrushes,

starting with finger brushes to electric brushes and then to innovative ultrasonic toothbrushes. Unlike electric brush models, a modern ultrasonic toothbrush also cleans the interdental spaces without any difficulties. It thoroughly removes bacteria and dental plaque, which can lead to gingivitis, dental calculus and resulting complications. Such deep cleaning is only possible due to innovative ultrasonic technology from Germany.

04

2.



More than

500.000
people

worldwide are already
enthusiastic users of
emmi®-pet.

ny
sny



OUR ULTRASONIC TECHNOLOGY

is powered by a piezoceramic chip located in the attachment, that moves with up to 96 million oscillations per minute.

emmi®-pet ultrasonic – good health starts with good teeth

emmi®-pet ultrasonic is dental care which appeals to cats and dogs. The patented hygiene system for teeth and mouth was especially developed to fulfil the needs of cats and dogs, and cleans teeth and gums silently and without any motions. Simply put a small amount of ultrasonic toothpaste on the brush attachment and hold it onto the teeth of your pet – emmi®-pet handles the actual cleaning on its own. Dental plaque and bacteria are removed gently without any scrubbing or humming, by 96 million ultrasonic air vibrations per minute. Via the ultrasonic vibrations contaminations are transported outwards and thus teeth and gums are cleansed thoroughly and gently – no chance for bad breath, dental calculus, gum diseases and resulting complications caused by bacteria in the mouth.

emmi® PET

The best dental care

starts when your dog is a puppy, in order to grow accustomed to dental cleansing as soon as possible.

The benefits of emmi[®]-pet ultrasonic



3.



PREVENTATIVE CARE TO THE FULLEST EXTENT

A regular use saves the additional visit to the vet for dog and owner: emmi®-pet ultrasonic removes dental plaque easily and prevents dental calculus as well as tooth and gum diseases*.

PATENTED ULTRASONIC

emmi®-pet ultrasonic is a patented hygiene system for teeth and mouth. In comparison to other electric toothbrushes (with about 40,000 mechanical motions) the piezo chip in the head of the brush cleans with up to 96 million ultrasonic air oscillation movements per minute.

WHITE TEETH

Via ultrasonic, the teeth stay shining and white.

SILENT AND WITHOUT ANY VIBRATIONS

emmi®-pet ultrasonic will not distress your pet.

HEALTHY TEETH, LIFE-LONG

If you care for teeth and gums regularly with emmi®-pet ultrasonic, the teeth will stay healthy the whole life long**.

FRESH BREATH

emmi®-pet ultrasonic removes germs and bacteria, which lead to bad breath. Meaning: bye-bye to the little stinkers!

* Dependent on how and how far a tooth is already affected.

** If the teeth are healthy and only the gums are affected, even a healthy tooth might become loose and fall out.

Helpful tips & tricks

around the dental health of your pet you will find on our website:

> www.emmi-pet.de

4



**Ultimately
suitable for
tooth and
fur care**



Ideal für
100% Ultraschall®
Made by EMMI Germany



EMMI®-PET TOOTH CARE ATTACHMENTS

The specially customized attachment for the thorough and gentle tooth care for cats and dogs is available in two sizes: the small attachment is suitable for cleaning incisors and little back teeth. The bigger attachment cleans canine teeth and bigger back teeth until they shine brightly. Cleaning the teeth one or two times a week with ultrasonic is completely sufficient for a strong and healthy set of teeth!

 **MADE IN GERMANY**

EMMI®-PET SKIN AND FUR CARE ATTACHMENTS

emmi®-pet has developed the Care ultrasonic attachment for gentle skin and fur hygiene. They are easily attached to the handle and moistened with water. Combined with the EICKEMEYER® ultrasonic gel you can ensure a positive effect on multiple symptoms like inflammation, tick bite or arthrosis by sonification with 84 million air vibrations per minute.

Shiny happy emmi®-pet customers and pets



SENSATIONAL RESULTS

"Just about one month we are using emmi®-pet now. Our Pluto got accustomed pretty well to brushing his teeth in the meantime. The results are sensational. I do not regret following the recommendation of a friend! The dental calculus of our dog has disappeared completely after only a few applications and even the gums have regenerated. Great compliments to the makers – simply great!"

Dietmar Renner



ASTONISHING

"I never would have thought this might work – but the dental calculus of my Labrador is gone. We fought a long time against it – bought special food etc. Unfortunately nothing seemed to work until I found emmi®-pet. We are using the equipment now for 5 months and are completely satisfied. To all dog owners – this really is a GOOD INVESTMENT for your dog! I would make this purchase again and again."

Anna Christina Stein



I DO NOT WANT TO MISS IT ANYMORE

"My daughter made emmi®-pet a present for me and my pug. I did not know the product at that time, but I do not want to miss it since the first treatment. The teeth of my pug are as healthy as never before. Even the vet was impressed by the results. I can only recommend emmi®-pet."

Jürgen Wolf



REALLY RECOMMENDABLE

"I bought the tooth brush about a year ago. In the beginning we brushed the teeth every day for three weeks. My Sally did not have any problems with it. I can do anything with her. After that we were able to scrape off the dental calculus. From this point on she was free of dental calculus. Now we are practicing the dental care one or two times a week. Do not be alarmed by the price. It really is worth it."

Jessica Egger





Toothpaste

emmi®-pet Basic Set
incl. charging station,
replacement brush &
toothpaste

Tooth care attachments

Fur and skin care
attachments

Emmi Ultrasonic GmbH



Emmi Ultrasonic is all about ultrasonic: The German-based company has specialized in development, manufacturing and worldwide distribution of innovative ultrasonic equipment and procedures. The patented emmi® ultrasonic products stand out by extremely miniaturized ultrasonic chips and innovative technologies that have been used for cleaning, preventative care and therapy in the oral and dermatological field since 2008. Emmi Ultrasonic is part of the EMAG AG Group, a medium-sized, family-owned company based in Moerfelden-Walldorf near Frankfurt am Main. All products made by EMAG AG are manufactured in the Rhine-Main area and are exemplary for the quality of the successful brand „Made in Germany“.



emmi[®] PET

Emmi Ultrasonic GmbH
Gerauer Strasse 34
D-64546 Moerfelden-Walldorf
Germany

Phone: +49 6105 406 794
Fax: +49 6105 406 781
info@emmi-ultrasonic.com
www.emmi-pet.com

 /emmipet
 @emmi_pet.com