

emmi®-skin Ultrasonic Crèmes & Foams

Ultrasonic Hyaluron crème reduces the formation of folds, tightens the skin, provides it with much needed moisture and induces the product of the body's own production of hyaluron. The HydraLift crème with apple stem cells serves as an anti-aging product. It slows the aging process, while regenerating and optimizing the skin's own production of collagen. The Ultrasonic Silicon crème counteracts the formation of folds by strengthening the connective tissues. Furthermore, it stimulates and promotes the regeneration of your own skin.



The Thermo-effect foam acts as heat therapy for joint and muscle pain, increasing the skin's vitality and promoting your circulation. It is worth noting that emmi®-skin crèmes will also work without emmi®-skin INTENSE ultrasonic unit. Of course, the emmi®-skin INTENSE will increase the positive effect of the crèmes.



Replacement Covers

You should be aware that the use of the ultrasonic cover by multiple people may not be appropriate for reasons of hygiene and it should be replaced after roughly three months of use. Only one person should use a given ultrasonic cover. You can purchase replacement covers at the same place where you purchased your emmi®-skin INTENSE unit.

Battery Charging Technology

The hand unit can be replaced on the charging station after use to re-charge the lithium battery or to be able to use the emmi®-skin INTENSE unit with a full battery charge every 8 days. If the period of usage is more than 6 minutes per day, the battery should be recharged for a brief period. When the battery is almost empty, the operational indicator (lowest green LED) will not light up. Instead, the red LED will light up, blink for 5 seconds and the unit will be completely turned off. When that happens, the unit must recharge for 24 hours so that it will work as expected.

Replacing the Battery & Protecting the Environment

The unit has a special, long-life, environmentally friendly lithium battery located in the hand unit. If the red LED continues to light up after long periods of use, despite recharging, the battery must be replaced. The battery can only be replaced by an emmi®- skin INTENSE lithium power-saver battery (original replacement part with an internal fuse). We will be happy to replace the battery for you, if you do not want to do this. If you want to replace the battery yourself, follow these instructions.

- 1 Loosen the screw at the bottom of the hand unit and pull the battery holder out with the floor plate.
- 2 When inserting the new battery, check the polarity. Match the battery's positive terminal with the floor plate holder.
- 3 Be sure that all of the seals sit properly when re-assembling the parts. Doing this ensure the IP54 class protection. You can purchase replacement batteries at the same place where you purchased your emmi®-skin INTENSE unit.

Note: any guarantee claims will only be applicable if the product has been treated properly and original parts have been used when replacing the battery.

Disposal of Batteries

Do not put the special battery in the unit in the household garbage. To dispose of the battery, please use the mutual return collection containers or other return centers that are available (such as any point of purchase for similar batteries). Please note the dangers indicated below.

Dangers



The emmi®-skin INTENSE unit includes a battery with a built-in fuse. Despite this, do not give the battery to children and comply with the notes about disposal above. If liquids leak from the battery despite the safety measures taken, do not allow them to contact the skin and dispose of the battery as described above. The connection with wall socket also requires safety measures. The charging station has been cast to insulate against electricity, however the charging station should never be placed in water. Make sure that the power cable for the charging station is free of all problems. If the housing, plug or power cable becomes damaged, the entire charging station should not be used anymore and will have to be replaced. The emmi®-skin INTENSE unit is only appropriate for the usage described above.

Guarantee

The guarantee will only apply to the 24 months after the date when the emmi®-skin INTENSE was purchased; the battery guarantee only applies for 6 months after the date of purchase. Defects in the form of material or manufacturing defects that appear during the guarantee period will be corrected free of charge. The correction may take the form of repairs or replacement. To the extent legally allowed, all other claims for compensation shall be excluded. The guarantee will be voided if the unit is opened by a third party or disassembled without the prior permission of Emmi Ultrasonic GmbH.

The guarantee will also be voided by intentional damage or incorrect usage. The guarantee period does not imply the provision of a guarantee. We retain the right to bill for any expenses that arise from unjustified complaints, such as errors made during installation or operation. The purchase receipt with the date of purchase must be kept to make guarantee claims. So that we can process your guarantee claim quickly, please the entire unit, including the ultrasonic cover and the charging station with sufficient postage directly to Emmi Ultrasonic GmbH. Shipments that are not free of charges will not be accepted. Thank you for your understanding.

Proper Usage

The products may only be used for the issues intended for in the catalogs and according to the technical descriptions, meaning cleaning and skin care purposes, and only in connection with the equipment and components recommended and allowed by Emmi Ultrasonic GmbH. ⚠ Any usage other than this or beyond this shall be considered as unintended. Emmi Ultrasonic GmbH will not assume any liability for improper usage, any damage arising and does not guarantee the problem-free function of the components. Compliance with the users' guide is a component of proper usage. The symbol on the name plate also refers to this. The smooth and safe operation of the products assumes proper transportation, storage, placement and assembly, as well as careful operation and maintenance. Compliance with repetitive checks in accordance with BGV A2 is mandatory. The unit must be setup so that the charging station's plug will be freely accessible at all times.

Improper Usage

Danger to life and limb as well as material damages may be caused by improper usage of the product. Improper usage is forbidden for that reason. Emmi Ultrasonic GmbH will not assume any liability for improper usage, any damage arising and does not guarantee the problem-free function of the components. Compliance with the following instructions is mandatory to exclude improper usage of the unit.

- The emmi®-skin INTENSE unit for skin cleaning using Original 100% emmi®- skin Ultrasonic may only be operated using the wall voltage that has been printed on the name plate (on the bottom of the charging station).
- Do not use combustible or aggressive chemicals, acids and similar agents for cleaning.
- If damage occurs during transportation do not connect the charging station to the wall.
- Only use the unit in internal rooms at room temperature.

emmi®-skin

INTENSE

Users' Guide

100% Ultrasonic for Young, Beautiful, Healthy Skin!



emmi ultrasonic

Emmi Ultrasonic GmbH • Gerauer Straße 34
64546 Mörfelden-Walldorf • Germany
Tel. +49 (0)6105 40 67 94 Fax +49 (0)6105 40 67 81
www.emmi-ultrasonic.de info@emmi-ultrasonic.de



EMAG AG
DIN EN ISO 9001:2008



RoHS CE
(S6Abs.3, ElektroG)

PATENTED, TESTED TECHNOLOGY
MADE IN GERMANY

**emmi®-skin INTENSE –
100% Ultrasonic for Young, Beautiful, Healthy Skin**

Congratulations! You have chosen a high-quality, innovative product. To always achieve the best possible results with your new emmi®-skin Ultrasonic Unit, please read the following users' guide carefully.

The emmi®-skin INTENSE Set includes:



1 emmi®-skin INTENSE hand unit, 1 Size L ultrasonic cover, 1 charging station, 1 bottle of Ultrasonic Crème (30 ml), 1 thermal heating sponge

Before Turning the Unit On

Plug the charging station into a wall socket. Place the hand unit on the charging station and charge it for roughly 24 hours before using it for the first time. If it is used daily for roughly 6 minutes, the battery charge will last roughly 8 days. Afterwards, re-charge the unit for 24 hours. The green LED will light up when the charging station connection has been made properly.

Before Initial Usage

Two additional green LEDs on the hand unit indicate the battery charging state. If your emmi®-skin INTENSE ultra sound unit has been charged properly, place the cover perpendicularly over the hand unit from above. When removing the ultrasonic cover, pull it off perpendicularly as well.

Activation and Deactivation

Activation:
Activate the unit by pressing the power button (for roughly 1 second). Optionally, you can deactivate the tone indicating a change in position by pressing the power button for roughly 9 seconds. The green operational indicator LED will light. The battery charging state will be indicated by two other blinking green LEDs. The unit will also vibrate as an additional functional indicator for the user.

Deactivation, Automated or Manual

Deactivate the unit by pressing the power button (for roughly 1 second). The unit will automatically turn off after 6 minutes (Power Saver mode for shipping). The unit also has an electronic system for detecting the cover. If the ultrasonic cover has not been put in place, is incorrect or defective, the unit will automatically deactivate and the lower green LED will turn off.

LED Indicators

1 Green POWER LED
Press the power button for roughly 1 second (activation protection) to turn the unit on. The lowest green LED will light (operational indicator).

2 Green Fully Charged LED
This blinking green LED indicates that the battery has been fully charged.

3 Green Partially Charged LED
This blinking green LED indicates that the battery has been charged to the 30% level.

4 Blue Charging LED
If power is being properly supplied to the hand unit at 100-240V, this LED will be blue. Once the unit has been completely charged, the LED will be darker (charging state check).

5 CHECK/FAILURE Red LED
If the battery is empty or there is a failure in the unit, this LED will be red. When the unit turns off, the red LED will blink for 5 seconds and then turn off.



Using the emmi®-skin INTENSE with 100% Original emmi®-skin Ultrasonic

When the unit has been turned on, you will feel it vibrate and hear a soft indicator tone. Because ultrasound cannot be seen, heard or felt, a small motor has been built into the hand unit so that the user can feel and hear when the unit has been activated and is ready for use. A soft tone will be generated every 6 seconds over the ideal usage period as well as when its position is changed. The motor will also stop briefly so that you receive tactile feedback. The tone can be turned off, as described above.

Usage

We recommend that you use the emmi®-skin INTENSE 100% Original Ultrasonic every 2 days for 6 minutes to clean the primary skin areas.

1 For reasons of hygiene, we recommend cleaning the ultrasonic cover before each use. For this purpose, we recommend using our disinfecting spray from the Emmi®-Beauty line, any other disinfecting spray or simple water.

2 Rub the unit over the primary skin areas to be treated, or combine the treatment with our emmi®-skin crème. The effect of the Original 100% emmi®-skin Ultrasonic unit is only guaranteed in combination with the use of emmi®-skin crème. You will find more information about the Ultrasonic Crème from emmi®-skin in the emmi®-skin Ultrasonic Crèmes & Foams section of this users' guide.

3 Turn the unit on and hold it against the location for at least 6 seconds without moving it. A longer application period will be harmless, because the emmi®-skin INTENSE uses a soft level of ultrasonic. We recommend an application period between 6 and 60 seconds for each primary area and problem. Longer application can lead to faster, better results.

4 Move to the next primary skin area and repeat Step 3.

5 Once the primary cleaning process has been completed, the emmi skin INTENSE unit will turn off automatically and. Clean the ultrasonic cover and the hand unit carefully with some water or a disinfecting spray. Afterwards, you can return the emmi®-skin INTENSE unit to the charging station with the cover on the unit.

What Our Customers Are Saying:

||||| A 92%
tighter skin feeling

||||| A 94%
improvement in skin appearance

||||| A 91%
increase in healthy appearance

* A 28- day test conducted daily with 24 test participants between 22 & 56 years of age. The emmi®-skin Hydralift, Hyaluron and Silicon Crèmes were used. The test participants indicated that their skin felt tighter (92%), that the appearance of their skin improved (94%) and that prior use of a crème in combination with the emmi®-skin INTENSE unit worked better (91%).

Effect

The emmi®-skin INTENSE unit has been developed by German engineers according to the most recent scientific discoveries. Unlike other traditional skin care devices, emmi®-skin INTENSE does not work by means of mechanical friction or abrasively like peeling. Instead, it uses the most modern ultrasonic technology, which does not abrade the skin. This technology has been used with wonderful success in many medical areas. Using ultrasound for skin care is absolutely harmless to people. The fact that using ultrasound at up to 50 watts is harmless to the human body has been proven scientifically. The emmi®-skin INTENSE uses 0.2 watts.

Dermatologists have tested the effects and safety of the technology. The patented piezo chip in the cover generates the ultrasonic waves and cause millions of micro-bubbles in the emmi®- Ultrasonic crème, which are channeled deep into layers of the skin. Human skin consists of 3 functional layers:

- the epidermis,
- the dermis,
- the hypodermis.

Conventional crèmes only affect the epidermis. In contrast, the emmi®-skin INTENSE Ultrasonic unit successfully channels the contents of emmi®-skin crèmes into the hypodermis and the leather of the skin due to its unique effect. Doing this provides the skin with effective support as it forms new cells and regenerates itself. The positive effect of using ultrasound also releases the body's own energies, improves the transportation mechanisms of the lymph nodes and improves circulation.

People are often told that water is required as a transfer medium. The reality is that the appropriate transfer medium depends on the frequency and performance rate of the ultrasonic waves. Medical applications are well aware of these issues and therefore use salves and gels with different consistencies as appropriate transfer media. We have developed a special ultrasonic crème for reducing the formation of folds and skin blemishes. To achieve the ideal results, please use emmi®- Crème together with emmi®-skin INTENSE.

